

BOFAP CARDBOARD BACTERIAL FILTER - MOUTHPIECE for Cosmed



Product Code	33910
Unit of sale	box of 200 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2370647
CND	R049001
EAN/UPC	8023279339109

Description

BOFAP BACTERIAL FILTERS - disposable

Internal x outer Ø: 28 x 30.5 mm

Compatible with: Cosmed

Minimum order: dispenser of 200 pcs.

BOFAP is a patented viral bacterial filter made of cardboard for spirometry examinations.

It combines the convenience of the mouthpiece with the usefulness of the filter.

It reduces the risk of cross-contamination among patients and helps to maintain equipment hygiene.

A bag of filter material is welded inside the mouthpiece to stop the same quantity of particulates of BV filters.

Laboratory tests have verified the ability to catch up to 99.99% of particulate matter as small as 0.2 microns in size.

BOFAP is fully recyclable and reduces plastic usage by 99%.

- reduced risk with lower costs than traditional plastic filters
- filtration efficiency comparable to plastic filters (above 99.99%)

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR

Made in Italy.