

# 3.5 MHz CONVEX PROBE for code 33965-6, 33877, 33882-3



Product Code	33953
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2975208
CND	Z11040182
EAN/UPC	6945121400107
UMDNS	15976

## Description

### 3.5 MHz CONVEX PROBE

- Scanning method: convex
- Frequency frequency (MHz): 3.5
- Probe band width QBit5 (MHz): 2-6.8
- Main application: abdominal, ob/gyn, urology