

3.0 MHz MICRO-CONVEX (BASIC CARDIAC) PROBE for 33960, 33877, 33882-3, 33965-6



Product Code	33964
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2975208
CND	Z11040182
EAN/UPC	6945121400367
UMDNS	15976
PARAF	903630022

Description

Probe for ecocolour Dopplers

- Scanning method: micro convex (BASIC CARDIAC)
- Main frequency: 3.0 MHz
- Probe band width QBit5 (MHz): 2-6.8
- Main application: CARDIOLOGY, PEDIATRICS