

2-5.0 MHz CONVEX PROBE for DP-30, DP-50, DP-50 EXP., DP-6600, Z5, Z50



Product Code	33982
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	170150
CND	Z1104018001
EAN/UPC	6944904001272

Description

3.5 MHz Convex probe.

- Mindray code: 35C50EA
- Scanning method: Convex
- Main freq.: 3.5
- Probe band width: 2.0 - 5.0

Main application:

- ABDOMINAL
- OB / GYN
- UROLOGY