

# 5-10 MHz LINEAR PROBE for DP-10, DP-20, DP-30, DP-50, DP-50 Expert



Product Code	33983
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	170156
CND	Z1104018001
EAN/UPC	6944904051567

## Description

Mindray Code: 75L53EA  
Scanning Method: Linear  
Main Frequency: 7.5 MHz  
Probe Bandwidth: 5.0 - 10.0 MHz  
Applications: Abdominal, Small Parts, MSK, Vascular, Pediatric, Intraoperative