

8-14 MHz LINEAR ARRAY PROBE for DP-30, DP50, DP50 Exp., Z5, Z50



Product Code	33990
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	317334
CND	Z1104018001
EAN/UPC	8023279339901

Description

10.0 MHz Linear array probe.

- Mindray code: 10L24EA
- Scanning method: Linear array
- Main freq.: 10.0
- Probe band width: 8.0 - 14

Main application:

- OB / GYN
- SMALL PARTS
- MSK
- VASCULAR
- PEDIATRIC

Not for DP-50PT