


NEW MINDRAY DP-50 EXPERT ULTRASOUND



non include
not included
non includes

Product Code	33992
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1947360
CND	Z11040104
EAN/UPC	6936415925796
GMDN	40761

Description

MINDRAY DP-50 EXPERT PORTABLE ULTRASOUND

With high resolution LCD monitor, 2 probe connectors, 500 GB hard disk, 4 USB ports, multiple software package, without probes.

Imaging mode: B, B/B, 4B, B/M, M

Image processing technology:

- Phase Shift Harmonic Imaging
- Tissue Harmonic Imaging
- iBeam™ Spatial Compound Imaging
- iClear™ Adaptive Speckle Suppression Imaging
- iTouch™ Auto Optimization
- Trapezoid imaging
- iZoom™ Automatically expand the image to full screen
- Slant scanning for linear probes (2D Steer)

Wideband multifrequency probes: continuously adjustable 2-14 MHz

Intelligent Workflow/Storage/Review:

- iStation™ Intelligent patient data management system
- 500 GB integrated hard disk for storage
- 4 USB ports
- Cine memory capacity (Max.) B mode: 10566 frames / M mode: 66.3 s

Supplied with GB, IT manual and CD manual GB, IT, FR, DE, RU, ES, PT, BR.

Technical Specifications

Display: 15" LCD
 Display language: GB, FR, IT, ES, PT, RU, PL, CZ, TR, FI, DK, NO, SE, DE
 Power: AC 100-240 V - 50/60 Hz
 Rechargeable lithium-ion battery - on request
 Size and weight: 14.7x36.1x35.8 cm (7.5"x16.3"x14.9"); 8.1 kg (17 lb)