

# 3.5 MHz CONVEX PROBE for ECO1/2/3/5/6



Product Code	33998
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2975214
CND	Z11040182
EAN/UPC	6945121407366
UMDNS	15976

## Description

Probe for Chison portable ultrasounds:

- Chison code: C3-A
- Scanning method: Convex
- Main Frequency: 3.5 MHz
- Frequency Range (MHz): 2,5 / 5,0
- Applications: ABDOMINAL, OB&GYN, UROLOGY, PEDIATRIC