


5.3 -10 MHz LINEAR PROBE for ECO1/2/3/5/6



Product Code	33999
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2975214
CND	Z11040182
EAN/UPC	6945121407397
UMDNS	15976

Description

Probe for Chison portable ultrasounds:

- Chison code: L7M-A
- Scanning method: linear (with 40 mm footprint, general linear probe)
- Main Frequency: 7.5 MHz
- Frequency Range (MHz): 5.3-10.0
- Applications: Small Parts, MSK, Vascular, Pediatric, Orthopaedic