

5-10 MHz LINEAR PROBE for DP-30, DP50, DP-50 Expert, Z5, Z50



Product Code	34004
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	170156
CND	Z1104018001
EAN/UPC	6944904037202

Description

7.5 MHz Linear probe.

- Mindray code: 75L38EA
- Scanning method: Linear
- Main freq.: 7.5
- Probe band width: 5.0 - 10.0

Main application:

- SMALL PARTS
- MSK
- VASCULAR