



|               |                   |
|---------------|-------------------|
| Product Code  | 34008             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8023279340082 |

## Description

CHISON ECO5 VET Colour Doppler 12" high resolution colour LED monitor, 8G high speed memory, 3 USB ports, 2 probe (not included) connectors, EASYVIEW™ archive system.

Imaging model: B, B/B, 4B, B/M, M, CFM, PW

Imaging processing:

- MCI (Multiple Compound Imaging) - trapezoidal
  - i-image intelligent image optimization - automatic PW trace
  - THI (Tissue Harmonic Imaging) - Chroma
  - SRA (Speckle Reduction Algorithm) - full screen display model
  - full digital beam-former - AIO (Automatic Image Optimization)
- Other specifications as ECO3EXP.

## Technical Specifications

- Dimension: 33.5 x 15.5 x h 35 cm (13.2"x6.1"xh 13.8")
- Weight: 6.5 kg
- Display: 12" LED, high-resolution colour 1,024 x 768
- Power:
  - AC 100-240 V, 50-60 Hz
  - Rechargeable lithium-ion battery

## Standard accessories

- 12" high resolution LED monitor - colour
- 8G memory card
- Hot key
- Pc printer connection
- 3 uSB port
- 1 ethernet port
- 1 video out port, 1 vGA out port
- 1 footswitch port
- Full screen display
- 2 hours rechargeable lithium-ion battery
- Power adaptor
- Standard software package
- Easy view image archive system
- THI
- PW (doppler function to detect blood flow speed)