

2 PIECES BELT - plastic lock - black



Product Code	34011
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1184395
CND	V0880
EAN/UPC	8023279340112
GTIN	08023279340112

Description

Immobilisation Belts for stretchers and spinal boards

Type: Belt-B
Colour: black
Size: 5 x 213 cm
Max load: 60 Kg
Buckle style: 2 swivel speed clip + plastic snap-lock
To be used with: spinal boards with pin