


SPINAL IMMOBILIZER - adult M



Product Code	34032
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2598639
CND	M03050201
EAN/UPC	6974442059138

Description

Spinal immobilizer - Adult M (150-175 cm height)

Designed to secure the patient, the surrounding is suitable for individuals with rib injuries. Strips, available in different colours, can be quickly fixed. It is equipped with nylon buckles and safety belts. It is recommended that when used in conjunction with a neck fixation, the patient should be in sitting position, before being transferred to a PE stretcher.

Supplied with forehead strap for head fixation and safety strap for chin fixation.
Max. load capacity: 175 kg.
Weight: 3 kg.
Unfolded size: 79x82 cm.
Folded size: 84x25x9 cm

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic