



Product Code	34040
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	5707480133857

## Description

Professional training manikin made to fulfil all the needs of Advanced Life Support (ALS) training. The focus of the manikin is on the algorithm: Airway Management, ECG simulation, defibrillation and the training of i.v. access technique.

Integrated Ambu Wireless Technology: the browser-based Ambu Manikin Management Module allows the wireless control, documentation and monitoring of the trainee's performance through the central interface, where all data and information come together: compression depth, correct hand position, ventilation volume, stomach inflation.

Customised scenario training: the intuitive user interface allows a quick and easy creation of custom scenarios and adapt them to different course concepts, in accordance to the demand of instructors and the current ERC / AHA guidelines. The CPR parameters and manual events are documented wirelessly and included into the evaluation and debriefing.

Individual trainee evaluation in group: more efficient training, each participant receives his personal results (up to four participants).

## Technical Specifications

Mouth to Mouth ventilation: NON  
 Mouth to nose ventilation: NON  
 Bag-valve-mask ventilation: YES  
 Use of Guedel airways: YES  
 Airway management: Et tube, laryngeal tube & mask and other SGAs: YES  
 Jaw-thrust manoeuvre: YES  
 Gastric insufflation: YES  
 Ambu hygiene system: NON  
 Chest compression: YES  
 Max compression depth: 7 cm  
 Max insufflation: 1.2 litres  
 Hand position display: YES  
 Infinitely adjustable chest stiffness: YES  
 Pulse: electrical  
 Placement of defibrillation electrodes: YES  
 Defibrillation: YES  
 I.V. access (need IV trainer 34041): optional  
 Wireless CPR software: YES  
 ECG display: YES  
 Carrying bag: YES  
 Weight: 13.3 Kg