



Product Code	34042
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	5707480043835

Description

Adult-sized training manikin that provide an optimal lifelike representation of the human anatomy. Mechanical instrument gives you instant feedback on ventilation volume, stomach inflation, chest compression depth in mm and wrong hand position.

Make intubation possible with ET tube, Guedel, Laryngeal Mask and combi-tube.

Ambu Airway man includes torso with hooks for arms, supplied in track suit jacket, training mat and carrying bag.

Size: 80 x 35 cm
Weight: 13 kg

Technical Specifications

Mouth to Mouth ventilation: NON
Mouth to nose ventilation: NON
Bag-valve-mask ventilation: YES
Use of Guedel airways: YES
Airway management: Et tube, laryngeal tube & mask and other SGAs: YES
Jaw-thrust manoeuvre: YES
Gastric insufflation: YES
Ambu hygiene system: NON
Chest compression: YES
Max compression depth: 7 cm
Max insufflation: 1.2 litres
Hand position display: YES
Infinitely adjustable chest stiffness: YES
Pulse: manual
Placement of defibrillation electrodes: YES
Defibrillation: NON
I.V. access (need IV trainer 34041): optional
Wireless CPR software: NON
ECG display: NON
Carrying bag: YES
Weight: 13 Kg