



Product Code	34045
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	5707480134939
PARAF	950834236

Description

Training manikin for education in school of all procedures for basic life support resuscitation. It can be used for realistic training of external chest compressions and the placement of AED pads. The unique hygienic system prevents cross contamination when demonstrating and training of the principals of ventilation through mouth and nose.

Guideline conform: it allows to perform CPR training following the latest international CPR guidelines.

Compact and portable: it has a low weight and compact packaging.

Easy AED training: due to the closed shaped torso, defibrillation pads can be placed to perform a realistic AED training.

Technical Specifications

Mouth to Mouth ventilation: YES
 Mouth to nose ventilation: YES
 Bag-valve-mask ventilation: NON
 Use of Guedel airways: NON
 Airway management: Et tube, laryngeal tube & mask and other SGAs: NON
 Jaw-thrust manoeuvre: NON
 Gastric insufflation: NON
 Ambu hygiene system: YES
 Insufflation bag (Ambu hygiene system): 20
 Manikin face mask (Ambu hygiene system): YES
 Chest compression: YES
 Max compression depth: 6 cm
 Max insufflation: 0.4 litres
 Hand position display: NON
 Infinitely adjustable chest stiffness: NON
 Pulse: NON
 Placement of defibrillation electrodes: YES
 Defibrillation: NON
 I.V. access (need IV trainer 34041): NON
 Wireless CPR software: NON
 ECG display: NON
 Carrying bag: NON (carton box)
 Weight: 0.8 Kg