



Product Code	34046
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	5707480114344
PARAF	950834352

Description

Inexpensive training manikin with all the well-known and sought-after Ambu look-listen-feel features.

The closed shape of the torso with its anatomical correct landmarks enables students to learn the complete algorithm of Basic Life support, including the use of Automatic External Defibrillators (AED) and the placement of the AED Pads.

Open Airway: the airway is opened only when the head is tilted - there are realistic airway resistance and chest rise.

Hygienic System: unique patented hygienic solution as the Ambu Man with disposable head bags and washable face pieces.

Compact Bag: the head of the Ambu Uniman+ can be placed inside the manikin. Everything is packed in a small bag.

Technical Specifications

- Mouth to Mouth ventilation: YES
- Mouth to nose ventilation: YES
- Bag-valve-mask ventilation: YES
- Use of Guedel airways: NON
- Airway management: Et tube, laryngeal tube & mask and other SGAs: NON
- Jaw-thrust manoeuvre: NON
- Gastric insufflation: NON
- Ambu hygiene system: YES
- Insufflation bag (Ambu hygiene system): 25
- Manikin face mask (Ambu hygiene system): YES
- Chest compression: YES
- Max compression depth: 7.5 cm
- Max insufflation: 1.8 litres
- Hand position display: YES
- Infinitely adjustable chest stiffness: NON
- Pulse: NON
- Placement of defibrillation electrodes: YES
- Defibrillation: NON
- I.V. access (need IV trainer 34041): NON
- Wireless CPR software: NON
- ECG display: NON
- Carrying bag: YES
- Weight: 4.5 Kg