



Product Code	34047
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	5707480132621

Description

Cost efficient and reliable training manikin designed to fulfill all requirements for the BLS algorithm training.

The closed shaped torso with its anatomical correct landmarks allows trainees to learn and train the use of Automatic External Defibrillators (AED) and the placement of pads. It allows an optimal training of chest compression and ventilation with realistic resistance and expiration through mouth and nose.

The unique hygienic system prevents cross contamination. The correct compression depth and required ventilation volume can be monitored by the feedback indicator on the manikin.

Ventilation as in real life: the trainees will feel the realistic airway resistance and the chest will start rising when the head is correctly tilted and air blown into the manikin.

Size: 70 x 37 cm.

Technical Specifications

Mouth to Mouth ventilation: YES
 Mouth to nose ventilation: YES
 Bag-valve-mask ventilation: YES
 Use of Guedel airways: NON
 Airway management: Et tube, laryngeal tube & mask and other SGAs: NON
 Jaw-thrust manoeuvre: YES
 Gastric insufflation: NON
 Ambu hygiene system: YES
 Insufflation bag (Ambu hygiene system): 25
 Manikin face mask (Ambu hygiene system): YES
 Chest compression: YES
 Max compression depth: 7.5 cm
 Max insufflation: 1.2 litres
 Hand position display: YES
 Infinitely adjustable chest stiffness: NON
 Pulse: NON
 Placement of defibrillation electrodes: YES
 Defibrillation: NON
 I.V. access (need IV trainer 34041): NON
 Wireless CPR software: NON
 ECG display: NON
 Carrying bag: YES
 Weight: 6.0 Kg