



Product Code	34048
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	5707480133864

## Description

Adult-sized training manikin that provide an optimal lifelike representation of the human anatomy. Mechanical instrument gives you instant feedback on ventilation volume, stomach inflation, chest compression depth in mm and wrong hand position.

Make intubation possible with ET tube, Guedel, Laryngeal Mask and combi-tube.

PC connection of Ambu CPR Software Kit, Ambu CPR Printer or USB Interface.

Ambu Airway man includes torso with hooks for arms, supplied in track suit jacket, training mat and carrying bag.

Size: 80 x 35 cm  
Weight: 13 kg

## Technical Specifications

Mouth to Mouth ventilation: NON  
Mouth to nose ventilation: NON  
Bag-valve-mask ventilation: YES  
Use of Guedel airways: YES  
Airway management: Et tube, laryngeal tube & mask and other SGAs: YES  
Jaw-thrust manoeuvre: YES  
Gastric insufflation: YES  
Ambu hygiene system: NON  
Chest compression: YES  
Max compression depth: 7 cm  
Max insufflation: 1.2 litres  
Hand position display: YES  
Infinitely adjustable chest stiffness: YES  
Pulse: manual  
Placement of defibrillation electrodes: YES  
Defibrillation: NON  
I.V. access (need IV trainer 34041): optional  
Wireless CPR software: YES  
ECG display: NON  
Carrying bag: YES  
Weight: 13 Kg