



| | |
|---------------|----------------|
| Product Code | 34057 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2233133 |
| CND | V0899 |
| EAN/UPC | 8023279340570 |

Description

MULTIPOSITION AUTO STRETCHER - wheelchair

Unique design multifunctional automatic stretcher, easily adjustable at different positions for emergency treatment, also in difficult situations. Front and back fixation to lock the stretcher into the ambulance at a desired position. The upper part can become a wheelchair. The stretcher can be adjusted to 6 levels, fold back and forth according to practical needs. When bent forward, the back of the stretcher is supported by fixed pneumatic springs with strapless controlling up to 60 degrees angle. Mattress, belts, side rails, I.V. stand and fixing device are included.

Technical Specifications

Maximum size (LxWxH): 198x56x89 cm
 Minimum size (LxWxH): 198x56x33 cm
 Packing size: 199x58x35 cm
 Load: 159 kg
 Net weight, Gross weight: 40 kg, 46 kg
 Maximum angle of leaning forwards backwards: 0-90 degrees
 The height scope of touching down (4 grades): 93 cm - 82 cm - 70 cm - 35 cm