



| | |
|---------------|----------------|
| Product Code | 34058 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| RDM (NSIS) | 2327931 |
| CND | V9099 |
| EAN/UPC | 8023279340587 |
| PARAF | 977813928 |
| GMDN | 41630 |

Description

EVACUATION CHAIR - orange/chrome

Designed and built for emergency patient removal from confined spaces.
Folded it can be stored in most ambulances.
Manufactured from chromed steel with brakes on the front wheels.
Seat cover made of vinyl coated nylon, easily washable.
Folds in one movement, equipped with two quick-release patient restraints.

Made in Italy.

Technical Specifications

Height (open/closed): 98/78 cm
Width: 41 cm
Depth with extended handles: open 80 cm
Depth when folded handles: closed 31 cm
Weight: 12 kg
Load: 116 kg