

# WHEELCHAIR STRETCHER - 2 wheels



Product Code	34060
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2299193
CND	V0899
EAN/UPC	8034028012693
PARAF	903163285

## Description

Specially designed for confined areas such as narrow hallways or multiple landing stairways; anywhere access is seriously limited.

Made in light aluminium with waterproof orange sheet.

Two rear stationary wheels Ø 12 cm, two back handles and four front handles, 2 adjustable 30 cm, facilitate transport on most floor types.

2 quick lock safety belts immobilize patient during transport. Locking safety latch maintains chair position during use.

## Technical Specifications

Folded size: 91 x 52 x 19.5 cm

Open size: 91 x 52 x 71 cm

Wheelchair size:

- back height: 91 cm - width: 52 cm

- seat height: 48 cm - width: 52 cm

Loading: 150 Kg

Weight: 8.5 Kg