



Product Code	34063
Unit of sale	1 kit
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
EAN/UPC	8023279340631
GTIN	08023279340631

## Description

### DOCTOR'S KIT 1

#### STANDARD ACCESSORIES:

- Utility bag - blue
- Aeon A200 non contact infrared thermometer (GB, FR, IT, ES)
- Oxy-3 finger oximeter (GB, FR, IT, ES, PT)
- Kobe palm aneroid sphygmomanometer
- Trad stethoscope - Y black
- GIMA LED penlight
- Fast tourniquet