

FULL AUTOMATIC MULTIPOSITION STRETCHER



| | |
|---------------|----------------|
| Product Code | 34072 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 2233137 |
| CND | V0899 |
| EAN/UPC | 8023279340723 |
| PARAF | 903163323 |

Description

FULL AUTOMATIC MULTIPOSITION STRETCHER

Full automatic multiposition stretcher. Easily adjustable to 6 levels for any emergency situation. Backwards and forwards reclining, back and foot rest adjustable to different positions. Can be loaded into ambulance by one operator only and locked at desired position. Structure in high-strength aluminium alloy. Soft foam surface mattress, covered by special seamless welded PVC. Protective block side rust proof and easy to clean. Mattress, belt, side rails, IV stand and fixing device are included.

Technical Specifications

Maximum size (LxWxH): 198x55x93 cm
Minimum size (LxWxH): 198x55x30 cm
Packing size: 220x60x35 cm
Load: 159 kg
Net weight, Gross weight: 39 kg, 45 kg