



Product Code	34074
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2299344
CND	V0899
EAN/UPC	8023279340747
PARAF	903163335

## Description

<div>WHEELCHAIR STRETCHER - 4 wheels</div></div>

<div>Specially designed for confined areas such as narrow hallways or multiple landing stairways; anywhere access is seriously limited.</div><div>Made of light aluminium with waterproof orange sheet. Equipped with 4 wheels, 1 with brake two back handles and four front handles, 2 adjustable 30 cm, facilitate transport on most floor types. </div><div>3 quick lock safety belts immobilize patient during transport.</div><div>Locking safety latch maintains chair position during use.</div>

## Technical Specifications

<div>Backrest size: height - width: 91 cm - 52 cm</div><div>Seat size: height - width: 48 cm - 52 cm</div><div>Folded size: 91x52x19.5 cm</div><div>Open size: 91x52x71 cm</div><div>Max. load: 150 kg</div><div>Weight: 8.5 kg</div>