



Product Code	34077
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2210086
CND	V0899
EAN/UPC	8023279340778
PARAF	903163374

Description

WHEEL STRETCHER - foldable in 2

Wheels, outriggers and adjustable backrest. High strength aluminium poles and vinyl coated, waterproof nylon canvas. Light weight, small sized, easy to carry.

Technical Specifications

- Open size: 188 x 53 x 21 cm
- Folded size: 92 x 53 x 12 cm
- Weight: 8 kg
- Load: 150 kg