



Product Code	34085
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2534386
CND	T030399
EAN/UPC	8023279340853
UMDNS	10-414
PARAF	951564208
GTIN	08023279340853
GMDN	10415

Description

Suitable for protection against cold and heat. Silver side in keeps the warmth of the body, gold side out reflects sunbeams.

Radiotransparent
Size: 160x210 cm
Weight: 63 g

Instructions: GB, FR, IT, ES, PT, DE, HU, RO, SE, FI, GR, HR, Arabic