

# PUMP for Vacuum Mattress - aluminium



Product Code	34095
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8056389233821
---------	---------------

## Description

Aluminium pump for vacuum mattress

Capacity: 0.75 l.

Rotating handle.

Supplied with hose and fitting.

Made in Italy.