

GIMA STRETCHER 2



Product Code	34101
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	42689
CND	V0899
EAN/UPC	8023279341010
PARAF	908734837

Description

GIMA STRETCHER 2 - foldable in 2

Emergency stretcher is sturdy, but lightweight in construction, with aluminium frame and two solid iron cross-bars for heavy use. Reinforced nylon canvas, fluid repellent and non-tear.

Technical Specifications

- Open size: 203 x 50 x 14 cm
- Folded size: 96 x 50 x 17 cm
- Weight: 6 kg
- Load: 120 kg