

FIRST AID AIR CUSHION SPLINTS - set of 4 pcs.



Product Code	34106
Unit of sale	1 set
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	72803
CND	Z12069080
EAN/UPC	8023279341065
PARAF	903163499

Description

FIRST AID AIR CUSHION SPLINTS - set of 4 pieces

Indispensable instrument for first aid in case of accidents, cutting and crushing injuries. Prevent surface bleeding without stopping the circulation and secure painless transportation of the injured part.
The bandage is a transparent cover of PVC easy to clean and reusable, resistant against water and shifting weather conditions.
X-rays can penetrate and medical diagnosis can be made without removing the bandage.

The long tube (14.5 cm) with the valve make the inflation of the bandage easy, even by the injured part himself.

Set of 4 different pieces: full leg, foot/ ankle, full arm, hand/wrist splints in a 4 compartment bag, with inflating pump.

Made in Italy.