



Product Code	34108
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2598630
CND	M030599
EAN/UPC	6974442051026

Description

HEAD IMMOBILIZERS

FERMO 3 - child

Colour: blue

Unfolded size (cm): 26x34x14

Weight (kg): 0.5

Material: high density waterproof plastic

Adult and child models of head immobilizers, different in shape, material and colour. All are composed of two head blocks and two straps to prevent any head movement of the injured person during transport. Both sides have an ear hole. They can be used with scooping and PE stretchers.

Waterproof plastic coating makes it easy to clean and prevents bacterial growth

Multilingual manual: GB, FR, IT, ES, PT, DE, NL,DK, SE, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.