



Product Code	34110
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2992424
CND	M9001
EAN/UPC	8023279341102
PARAF	903117366
GTIN	08023279341102
GMDN	60790

## Description

INSTANT ICE - TNT 14x18 cm bag (non woven fabric)  
minimum order: pack of 25 bags (38x25xh 17 cm)

Icy-pack "instant ice" is recommended for cryotherapy.

Excellent in case of traumas to soothe the pain and to prevent haematomas, swellings, etc.

Easy and quick to use anywhere and anytime.

Press the bag to break the inner bubble and shake to mix contents. The blend you obtain drops temperature of the pack down. It reaches -4°C upon activation, it reaches 0°C in 20 minutes (at room temperature of 20°C).

Weight: 250 g.

Instant ice with multilingual bag: GB, FR, IT, ES.