



Product Code	34114
Unit of sale	box of 30 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2883469
CND	M9099
EAN/UPC	8683229609807
PARAF	913758797

Description

HOT-COLD THERMO GEL - box of 30 pcs.

Special gel keeps temperature (hot or cold) for long time.

Size: 14x28 cm.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.

Individually packed in single box.