

BODY BAG - vinyl - white - load 150kg



Product Code	34118
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279341188

Description

BODY BAG - vinyl - white
Made of vinyl with one "U" zipper and 4 handles.

Size: 210x90 cm; 0.02 cm thickness
Max loading: 150 kg

Multilingual manual: GB, FR, IT, ES, PT, DE, PL, DK, SE, FI, NO, NL, CZ, SK, HU, HR, SI, RO, BG, GR, LT, LV, EE, Arabic