

SPLINTS - size 15 x 495 mm



Product Code	34120
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2461414
CND	M030505
EAN/UPC	8023279341201
PARAF	903163537

Description

ALUMINIUM SPLINTS soft lining - 15 x 495 mm

Soft lining for comfort.

Aluminium thickness 1.0 mm.

Multilingual box and instructions: GB, FR, IT, ES, DE, HR.

Made in Italy.