



Product Code	34121
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2461414
CND	M030505
EAN/UPC	8023279341218
PARAF	903163552

Description

ALUMINIUM SPLINTS soft lining - 20x495 mm

Soft lining for comfort.
Aluminium thickness 1.0 mm.
Multilingual box and instructions: GB, FR, IT, ES, DE, HR.

Made in Italy.