

# SPLINTS - size 25 x 495 mm



Product Code	34122
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2461414
CND	M030505
EAN/UPC	8023279341225
PARAF	903163564

## Description

ALUMINIUM SPLINTS soft lining - 25x495 mm

Soft lining for comfort.  
Aluminium thickness 1.0 mm.  
Multilingual box and instructions: GB, FR, IT, ES, DE, HR.

Made in Italy.