



Product Code	34135
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No
RDM (NSIS)	138553
EAN/UPC	8034028018039

Description

Special kit designed for first aid in sporting clubs (tennis clubs, swimming pools, gymnasium,...).
Soft nylon bag has 2 side packets and zip closure.

Technical Specifications

- Size: 470 x 240 x 190 mm
- Weight: 6.5 kg

Standard accessories

- 1 Isothermic blanket Gold/Silver 160 x 210 cm
- 1 Disinfectant liquid 250 ml
- 2 Adhesive plasters 5 m x 2.5 cm
- 2 Pairs of sterile gloves
- 2 Vento mouth to mouth resuscitators
- 6 Sterile gauze pads 18 x 40 cm
- 2 Instant ice sprays 150 ml
- 3 Triangular TNT bandages 96 x 96 x 136 cm
- 2 Elastic bandages 6 cm
- 2 Elastic bandages 8 cm
- 1 Cotton wool pack
- 1 Pair of LISTER scissors 14.5 cm
- 2 Compressive pad bandages 80 x 100 mm
- 2 Compressive pad bandages 100 x 120 mm
- 1 Pack of 20 adhesive plasters 2 x 7 cm
- 8 Adhesive plasters 10 x 6 cm
- 1 Sterile TNT bandage 60 x 40 cm
- 1 Sterile TNT bandage 60 x 80 cm
- 6 Sterile TNT pads 10 x 10 cm
- 5 Instant ice packs
- 1 Hydrogen peroxide solution
- 1 Elastic net bandage
- 1 Haemostatic lacet
- 1 Sterile tweezers
- 5 Ammonia towelettes
- 1 Sterile saline solution bottle 250 ml
- 1 Gel hands disinfectant 100 ml
- Guidance leaflet 5 languages