



Product Code	34137
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No
RDM (NSIS)	136260
EAN/UPC	8034028011146

Description

SOFT BAG KIT DIN 13164-2014

Contents conform to DIN norms 13164.

Size: 260x172x82 mm.

Weight: 0.7 kg.

Standard accessories

1 adhesive plaster roll 5 m x 2.5 cm
 4 disposable vinyl gloves
 2 gauze bandages 6 cm x 4 m
 3 gauze bandages 8 cm x 4 m
 1 pair of lister scissor 14.5 cm
 1 isothermic blanket gold/silver 160x210
 1 sterile compressive pad 60x80 mm
 2 sterile compressive pads 80x100 mm
 1 sterile compressive pads 100x120 mm
 1 sterile TNT burn sheets 60x40 cm
 1 sterile bandage TNT 60x80 cm
 3 sterile gauzes 10x10 cm
 2 alcohol 70% wipes
 1 pack 14 assorted patches
 2 triangular bandages TNT 96x96x136 cm
 Instructions (GB, FR, IT, ES, DE, TR)