

# EMPTY METAL CABINET 370x300x140 mm



|               |                   |
|---------------|-------------------|
| Product Code  | 34141             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|        |              |
|--------|--------------|
| UK-REP | Not required |
| CH-REP | Not required |

|         |               |
|---------|---------------|
| EAN/UPC | 8023279341416 |
| PARAF   | 903163778     |

## Description

Empty metal cabinet

- Size: 370 x 300 x 140 mm