



Product Code	34143
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

RDM (NSIS)	141416
EAN/UPC	8056389239892

Description

Metal cabinet supplied with 3 internal metal shelves, removable bottle floor support, 2 small shelves on the door, lock and key.

Size: 402 x 582 x 202 mm
Weight: 12.9 kg

Standard accessories

- 3 Bags liquid soap
- 3 Bags antiseptic wipes
- 1 Bottle antiseptic solution 250 ml
- 3 Bottle antiseptic solution 500 ml
- 1 Oxygenated water bottle
- 2 Ammonia bag
- 5 Pair protective gloves, sterile
- 1 Box of 100 assorted plasters Plastosan
- 1 Box of 100 assorted 7 x 2 cm Plastosan
- 2 Adhesive plaster roll 2.5 cm x 5 m
- 1 Elastic bandage with clip 4.5 cm x 7 m
- 1 Pack of 8 assorted gauze bandages
- 1 Tubular net gauze
- 10 Pack of sterile compress 10 x 10 cm
- 6 Pack of sterile compress 18 x 40 cm
- 3 Pack of cotton wool
- 1 Mask with face shield
- 1 Pair of scissors cm 14.5
- 2 Sterile forceps
- 3 Haemostatic lacet
- 1 Sphygmomanometer + Stethoscope
- 2 Sterile blanket 40 x 60 cm
- 2 Pack of instant ice, single use
- 1 Digital thermometer
- 5 Trash bag
- 3 Irrigating solution 500 ml
- 4 Triangular blanket 96 x 96 x 136 cm
- 1 Isothermic blanket gold/silver 160 x 210cm
- 3 Sterile burn gel sachet 3.5 g

- 1 Compressive Emocontrol bandage
- 1 Pinocchio helicoidal open mouth Vento