

INTUBATION KIT with 3 standard laryngoscope blades



Product Code	34152
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP No

CH-REP No

EAN/UPC 8023279341522

PARAF 903163994

Description

INTUBATION KIT - complete

Compact nylon bag suitable for various uses: 2 internal zip pockets, 9 elastic retainers for instruments, made in water-resistant anti-scratch 600D polyester, blue colour.

Size: 32 x 22.5 x h 7.5 cm.

Standard accessories

- 3 Mc Intosh Laryngoscope s/s blades (1, 2, 3)
- 1 Laryngoscope metal handle
- 1 Magill forceps
- 1 Tongue holding forceps
- 3 Guedel cannula
- 1 Mouth opener
- 1 Scissors cm 14.5