

VENTURI TYPE VARIABLE CONCENTRATION MASK - adult

Product Code	34170
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2364599
CND	R03010202
EAN/UPC	8033003260425
GMDN	35175



Description

Disposable adult mask with Venturi system, with corrugated tubing 15 cm, extension tubing 200 cm, cup and seven adapters for variable oxygen concentration.

Latex free.

Technical Specifications

Colour code according to EN 13544-3: 2001

• Blue

Regulators: 22%

Liters / minute to reach the desired oxygen concentration: 2

• White

Regulators: 28%

Liters / minute to reach the desired oxygen concentration: 4

• Orange

Regulators: 31%

Liters / minute to reach the desired oxygen concentration: 6

• Yellow

Regulators: 35%

Liters / minute to reach the desired oxygen concentration: 8

• Red

Regulators: 40%

Liters / minute to reach the desired oxygen concentration: 8

• Fuchsia Pink

Regulators: 50%

Liters / minute to reach the desired oxygen concentration: 12

• Green

Regulators: 60%

Liters / minute to reach the desired oxygen concentration: 15