

"GIMA 7" EMERGENCY BAG



Product Code 34174

Unit of sale 1 pc

Minimum order 1

Type Kit of medical devices

Class

UK-REP No

CH-REP No

EAN/UPC 8023279341744

PARAF 977814033

Description

EMERGENCY KITS GIMA 7

Smart bag - large - red 600D polyester (27153)
Size: 65 x 35 x h 35 cm

Products included in kits:

- Emergency cover (34085)
- Universal scissors - 16.5 cm (20571)
- Facemask n. 2 - child/small (34222)
- Silicone resuscitator bag with mask n 4 - adult (34245)
- Teeth screw driver (34273)
- Tongue holding forceps (34274)
- Oxygen tube - 120 cm (34276)
- Mixed Guedel airways - box of 10 (34439)
- Oxygen therapy mask - adult (34166)
- 0.5 l oxygen cylinder (empty) with reducer UNI std (34501)
- Jotarap dual head stethoscope - black (32580)
- London sphygmomanometer (32725)
- Laryngoscope set 3 Mc-Intosh blades 2-3-4 (34302)
- Hand suction unit (28125)
- Delta pen light (25629)
- Fast tourniquet - green (25728)
- Magill forceps - 20 cm (34002)
- Gima digital thermometer °C (25560)
- Disinfectant wipes - tube 100 wipes (36625)
- Emocontrol elastic bandage (34782)
- Gauze bandage 3.5 m x 7 cm - box of 10 (34841)
- Sterile adhesive dressing 10 x 10 cm - box of 5 (34819)
- Oxy-4 fingertip oximeter (35092)
- Instant ice - bag 3 pcs. (34111)