

MEDIUM KIT - metal cabinet



Product Code	34197
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

RDM (NSIS)	142027
EAN/UPC	8034028010033
PARAF	910552811

Description

- Medium kit - metal cabinet
- Size: 370 x 300 x 140 mm
 - Weight: 3.8 kg

Standard accessories

- 3 Bags liquid soap
- 3 Bags antiseptic wipes
- 1 Bottle antiseptic solution 125 ml
- 2 Ammonia bag
- 2 Pair protective gloves, sterile
- 1 Box of 10 assorted plasters
- 1 Adhesive plaster roll 2.5 cm x 5 m
- 1 Gauze bandage 10 cm x 3.5 m
- 3 Pack of sterile compress 10 x 10 cm
- 3 Pack of sterile compress 18 x 40 cm
- 1 Pack of cotton wool
- 1 Pair of scissors 14.5 cm
- 1 Sterile forceps
- 1 Haemostatic lacet
- 1 Sterile blanket 40 x 60 cm
- 1 Pack of instant ice, single use
- 1 Trash bag
- 1 Irrigating solution 250 ml
- 1 Triangular blanket 96 x 96 x 136 cm