

AMBU MARK IV RESUSCITATOR - adult - with reservoir bag 1,500 ml



Product Code	34201
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	Kit of medical devices
------	------------------------

Class	
-------	--

UK-REP	Yes
--------	-----

EAN/UPC	5707480019434
---------	---------------

Description

AMBU® MARK IV RESUSCITATOR - adult - with reservoir

Designed in accordance with the original double wall principle, the new 100% latex-free Ambu bags are more robust and provide users with exceptional visual and tactile lung condition feedback during manual resuscitation, making operation safer and more secure. Built-in pressure limitation prevents compression if there is any resistance in the patient's airways.

Suitable for adults and children over 15 kg (3 years).

Supplied with face mask N°5, reservoir bag 1500 ml and single shutter valve.

Autoclavable at 134°C.
Maximum volume 1300 ml.