

AMBU MARK IV IN BAG RESUS. + accessories



Product Code	34217
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279342178
GTIN	08023279342178

Description

Designed in accordance with the original double wall principle, the new 100% latex-free Ambu bags are more robust and provide users with exceptional visual and tactile lung condition feedback during manual resuscitation, making operation safer and more secure.
Built-in pressure limitation prevents compression if there is any resistance in the patient's airways.

Suitable for adults and children over 15 kg (3 years).
Supplied with face mask N°5 and single shutter valve, helicoidal mouth-wedge, tongue holding forceps, 3 Guedel airways (small, medium and large) and oxygen tube, all fitted in a resistant nylon bag.
Autoclavable at 134°C.
Maximum volume 1300 ml.