

CPR MASK - pocket resuscitator



Product Code	34218
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2527128
CND	R03020201
EAN/UPC	4716770662404
UMDNS	15230
PARAF	903630376
GTIN	08023279342185

Description

CPR MASK - pocket resuscitator

The CPR pocket mask is a non-invasive, latex free mask with an integral non-rebreathing valve and efficient filter to prevent direct mouth contact with the patient's face when providing mouth-to-mouth resuscitation. It can be used on adult, child and infant. Supplied in a compact plastic box with a corrugated tube for an easier use.