



| | |
|---------------|----------------|
| Product Code | 34222 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|----------------|
| RDM (NSIS) | 2526739 |
| CND | R030199 |
| EAN/UPC | 8023279342222 |
| UMDNS | 12447 |
| PARAF | 903164236 |
| GTIN | 08023279342222 |
| GMDN | 36066 |

Description

GIMA AUTOCLAVABLE FACE MASKS

Made in silicone. Completely transparent, allowing a total visibility. They can be used with any resuscitator and lung ventilator.

Facemask N° 2 - child/small - latex free. Autoclavable up to 40 times.