

# AMBU SPUR II RESUSCITATOR - adult with connector for demand valve



Product Code	34234
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

GTIN 05707480043637

## Description

AMBU SPUR II ADULT - with face mask n° 5 and reservoir - connector for demand valve - 325023000

AMBU SPUR II RESUSCITATOR - single patient

Ambu SPUR II provides users with exceptional tactile and visual feedback during resuscitation. The bag is highly responsive, with minimal mechanical resistance. Furthermore the characteristic design provides optimum stroke volume with perfect recoil. Range includes infant, pediatric, and adult sizes.

All Ambu SPUR II resuscitators come in individual, resealable plastic carrying bags, complete with mask, reservoir bag, oxygen tube.

The bags are color-coded for fast identification of infant, pediatric, and adult sizes.

## Technical Specifications

Suitable for: adult > 30 kg  
Stroke volume 1 hand 600 ml / 2 hands 1,000 ml  
Resuscitator volume: 1,510 ml  
Size (length x diameter): 284x127 mm  
Bag reservoir volume: 2,600 ml